

STERILE

Klerwipe Sporidical Enhanced Peroxide

Klerwipe Sporidical Enhanced Peroxide pouch wipes are sterile, non-woven, low particulate, hydroentangled polyester/cellulose wipes impregnated with a blend of 1.5% hydrogen peroxide and deionised water (DI).

Liquid Specification

Colour	Colourless to slight yellowish	pH	2.01 - 2.30
Odour	Slight	Specific gravity @20°C	1.010 - 1.016
Clarity	Transparent liquid	Content of H ₂ O ₂	1.35 - 1.65%
Filtration	0.2 µm		

Wipe Specification

Physical Properties	Typical Values	Specification	Test Method
Material Composition	55% Cellulose/ 45% Polyester	55±5% Cellulose/45±5% Polyester hydro-entangled wipe	
Colour	White	White	Visual assessment
Basis weight (g/m ²)	68	68±3	Micro scales
Chemical Compatibility	Contains no chemical binders or optical brightener and is not affected by most common solvents & chemicals		

Extractables	Typical Values	Specification	Test Method
Particles (≥0.5µm) (biaxial shake)(x 10 ³ /cm ²)	16	≤45	IEST-RP-004.4
Fibres (≥100µm/cm ²)	200	≤410	IEST-RP-004.4
Non Volatile Residues: H ₂ O (g/m ²)	0.030	≤0.06	IEST-RP-004.4
2-Propanol (g/m ²)	0.034	≤0.04	
Extractable Ions (ppm) Sum of Ions (Na ⁺ , K ⁺ , Mg ²⁺ , Ca ²⁺)	69.8	≤500	IEST-RP-004.4

Packaging Information

Production & Packaging	Klerwipe Sporidical Enhanced Peroxide wipes are manufactured using Klercide Sporidical Enhanced Peroxide which is 0.2 µm filtered and sterile filled under Grade A laminar flow and double bagged in a Grade B cleanroom
Sterility	Preassembled dry goods are sterilised by gamma irradiation and Klerwipe Sporidical Enhanced Peroxide wipes are manufactured utilising sterile fill technology
Product Expiry	18 months from date of manufacture
Additional Testing	Supplementary data available to demonstrate 7 days chemical stability post package opening

Product Ordering Information

Scan or click the QR code or for general product ordering information visit www.ecolablifesciences.com/product



Data relates to the pre-impregnated dry wipes only. Please refer to the relevant Klercide Technical file for details of the impregnate.

TS064E

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

Please speak to your Ecolab account manager for further information or call 02920 854 390 (UK) +44 2920 854 395 (Export)

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences

EUROPE HEADQUARTERS

Richtistr. 7
8304 Wallisellen
Switzerland

www.ecolablifesciences.com

