



**BD**

Indispensable to  
human health

## **BD Sterile Pack Swabs** For Environmental Monitoring

Maximum sterility with double wrapped  
packaging, gamma-irradiation and  
performance-validation.

# BD Sterile Pack Swabs For Environmental Monitoring

Maximum sterility with double wrapped packaging, gamma-irradiation and performance-validation.



## DON'T TAKE CHANCES



that your environmental samples could be contaminated by a less than sterile swab. Wipe out uncertainty with BD Sterile Pack Swabs – the first pre-filled swab/rinse solution set with a Sterility Assurance Level (SAL) of  $10^{-6}$ , an industry standard! BD Sterile Pack Swabs are ideal for the stringent sterility requirements of surface sampling in clean rooms and isolators at pharmaceutical and medical device manufacturing facilities. And because these swabs are designed for hard-to-reach places, BD Sterile Pack Swabs easily fit into equipment recesses, nooks and crevices.

## PRODUCT FEATURES

- Packaged 50 to a box, with 10 single wrapped swabs in each double wrapped pouch
- Polypropylene tube is pre-filled with a 10 mL rinse solution
- Simple sampling and transport—easy-to-use swab is fixed to the screw cap
- Complies with the new ISO guidelines
- Each lot comes with a Certificate of Analysis
- ATP-free, Dacron™ swab
- Saves time, no need to create your own media and swab/solution set in house
- Convenient and cost-effective solution
- One year shelf-life, at room temperature storage

## PROCEDURAL STEPS



1. Open double wrapped exterior and peel open single wrapped pouch.



2. Twist cap and remove swab from tube.



3. Swab sampling area.



4. Twist cap onto tube and transport to the laboratory for testing.



**BD Diagnostic Systems**  
www.bd.com/industrial

## ORDERING INFORMATION

Cat. Number	Description	Qty/Pkg.
220518	BD Sterile Pack Swabs	200

### Local Offices

#### Austria

Am Concorde Park E1/7  
2320 Schwechat  
Tel: 43.1.70.63.66.00

#### Eastern Europe/ Middle East, Africa

Tullastrasse 8-12  
69126 Heidelberg  
Tel: 49.6221.305.184

#### Germany

Tullastrasse 8-12  
69126 Heidelberg  
Tel: 49.6221.305.0

#### Italy

Via Caldera 21  
20153 Milano  
Tel: 39.02.48.24.01

#### Spain

Camino de Valdeoliva s/n  
San Agustin de Guadalix  
28750 Madrid  
Tel: 34.91.848.81.00

#### Switzerland

Binningerstrasse 94  
4123 Allschwil  
Tel: 41.61.485.22.22

#### Belgium

POB 13, Denderstraat 24  
9320 Erembodegem  
Tel: 32.53.720.211

#### France

11 rue Aristide Bergès  
38800 Le Pont de Claix  
Tel: 33.4.76.68.36.36

#### Greece

5 Amfitheas Avenue  
17122 Nea Smyrni  
Athens  
Tel: 30.210.940.77.41

#### Netherlands

Postbus 757  
2400 AT Alphen aan  
den Rijn  
Tel: 31.20.582.94.16

#### Sweden

Instrumentvägen 19  
PO Box 32054  
12611 Stockholm  
Tel: 46.8.775.51.00

#### U.K.

21 Between Towns Road  
Cowley  
Oxford OX4 3LY  
Tel: 44.1865.781.677